

Grinder-Polisher

Tegrapol 35 / Tegraforce 5

Principle

Semi-automatic grinder and polisher for mounted samples with power head.

Manufacturer	Struers
Slice Diameter	300 mm
Rotational Speed Slice	40 - 600 rpm
Rotational Speed of Specimen Holder	50 - 150 rpm
Rotational Direction of Specimen Holder	Co-rotating / counter-rotating
Specimen Holder	Mounted specimen: 25 mm Mounted specimen: 40 mm Mounted specimen: 48 x 40 mm
Force Individual Specimen	10 N - 100 N
Force Specimen Holder	30 N - 400 N
Accessory	Tegradoser, automatic dosing
