

Cutting Grinder

Accutom 50

Principle

Cutting of small samples.

Manufacturer	Struers A/S
Revolution speed	300 - 5000 U/min
Rate of travel	0.005 – 3.00 mm/s
Positioning Range	Y-area: 105 mm X-area: 60 mm
Sample Sizes	Length of the sample max. 30 mm Width of the sample max. 225 mm
Circular Saws	Diameter: 75 to 152 mm max. thickness: 2 mm hole: ø 12,7 mm
Circulation cooling	Content: 4 Liter
