

Grinder-Polisher

Phönix 4000

Principle

Semi-automatic Grinder and Polisher for mounted samples with power head.

Manufacturer	Buehler Ltd.
Disc Diameter	300 mm
Rotational Speed	50 - 600 rpm
Rotational Speed of Specimen Holder	150 rpm
Rotational Direction of Specimen Holder	Clockwise/ counter-clockwise/ differential run
Specimen Holder	Mounted specimen: 25 mm Mounted specimen: 40 mm
Force Individual Specimens	- 100 N
