

Grinder-Polisher

Rotopol 31 / Rotoforce 4

Principle

Semi-automatic Grinder and Polisher for mounted samples with power head.

Manufacturer	Struers
Disc Diameter	300 mm
Rotational Speed of Disc	30 - 600 rpm
Cooling of Disc	Yes
Rotational Direction of Specimen Holder	150 rpm
Rotational Direction of Specimen Holder	Clockwise/ counter-clockwise/ differential run
Specimen Holder	Mounted specimen: 25 mm Mounted specimen: 40 mm Mounted specimen: 48 x 40 mm
Force Individual Specimens	5 N - 100 N
Force Specimen Holder	30 N - 400 N
